BPL_STEM_AIR_Perfusion - test

Here we show simulations of stem cell cultivation in an aerated hollow fiber reactor. The reactor volume is kept constant and cells recycled, thus the setup is similar to perfusion cultivation.

The model combines rudimentary cell growth and metabolism combined with times series data of of the metabolic rates: qNmax, qLc, and qO2, marked with red in the comprehensive plot.

Ref Greuel et al: "Online measurement of oxygen enables continuous noninvasive evaluation of human-induced pluripotent stem cell (hiPSC) culture in a perfused 3D hollow-fober bioreactor", Biotech. Bioeng., 2019.


Note that both disp() and describe() takes values from the last simulation and the command process diagram() brings up the main configuration

Brief information about a command by help(), eg help(simu)
Key system information is listed with the command system_info()

```
In [2]: plt.rcParams['figure.figsize'] = [30/2.54, 24/2.54]
In [3]: #process_diagram()

In [4]: # Process parameters to mimic stem cell example
    # - study impact of different values maintenatnce qm


newplot(plotType='Basic')
    for value in [1.0e-6, 10.0e-6, 20e-6]:
        par(qm=value)
        simu(350)
```



```
In [5]: # Process parameters to mimic stem cell example
    # . comprehensvie plot including metabolic rates

par(qm=1e-6)

newplot(plotType='Comprehensive')
simu(350)
```



Appendix

```
In [14]: system_info()

System information
   -OS: Linux
   -Python: 3.10.14
   -Scipy: not installed in the notebook
   -PyFMI: 2.14.0
   -FMU by: OpenModelica Compiler OpenModelica 1.23.1
   -FMI: 2.0
   -Type: FMUMModelME2
   -Name: BPL_STEM.Reactor
   -Generated: 2024-09-06T19:04:53Z
   -MSL: 3.2.3
   -Description: BPL - not used
   -Interaction: FMU-explore version 1.0.0
In []:
```