

✓ BPL_STEM_AIR_Perfusion script with PyFMI

The key library PyFMI is installed.

After the installation a small application BPL_STEM_AIR_Pefusion is loaded and run. You can continue with this example if you like.

```
!lsb_release -a # Actual VM Ubuntu version used by Google
```

```
↳ No LSB modules are available.
   Distributor ID: Ubuntu
   Description:    Ubuntu 22.04.4 LTS
   Release:        22.04
   Codename:       jammy
```

```
%env PYTHONPATH=
```

```
↳ env: PYTHONPATH=
```

```
!python --version
```

```
↳ Python 3.11.11
```

```
!wget https://repo.anaconda.com/miniconda/Miniconda3-py311_24.11.1-0-Linux-x86_64.sh
```

```
!chmod +x Miniconda3-py311_24.11.1-0-Linux-x86_64.sh
```

```
!bash ./Miniconda3-py311_24.11.1-0-Linux-x86_64.sh -b -f -p /usr/local
```

```
import sys
```

```
sys.path.append('/usr/local/lib/python3.11/site-packages/')
```

```
↳ --2025-03-25 16:18:26-- https://repo.anaconda.com/miniconda/Miniconda3-py311_24.11.1-0-Linux-x86_64.sh
   Resolving repo.anaconda.com (repo.anaconda.com)... 104.16.32.241, 104.16.191.158, 2606:4700::6810:bf9e, ...
   Connecting to repo.anaconda.com (repo.anaconda.com)|104.16.32.241|:443... connected.
   HTTP request sent, awaiting response... 200 OK
   Length: 145900576 (139M) [application/octet-stream]
   Saving to: 'Miniconda3-py311_24.11.1-0-Linux-x86_64.sh'
```

```
Miniconda3-py311_24 100%[=====] 139.14M 70.4MB/s in 2.0s
```

```
2025-03-25 16:18:28 (70.4 MB/s) - 'Miniconda3-py311_24.11.1-0-Linux-x86_64.sh' saved [145900576/145900576]
```

```
PREFIX=/usr/local
```

```
Unpacking payload ...
```

```
Installing base environment...
```

```
Preparing transaction: ...working... done
```

```
Executing transaction: ...working... done
```

```
installation finished.
```

```
!conda update -n base -c defaults conda --yes
```

```
↳ Channels:
   - defaults
   Platform: linux-64
   Collecting package metadata (repodata.json): done
   Solving environment: done
```

```
## Package Plan ##
```

```
environment location: /usr/local
```

```
added / updated specs:
  - conda
```

The following packages will be downloaded:

package	build	
ca-certificates-2025.2.25	h06a4308_0	129 KB
certifi-2025.1.31	py311h06a4308_0	163 KB
openssl-3.0.16	h5eee18b_0	5.2 MB
Total:		5.5 MB

The following packages will be UPDATED:

```

ca-certificates      2024.11.26-h06a4308_0 --> 2025.2.25-h06a4308_0
certifi              2024.8.30-py311h06a4308_0 --> 2025.1.31-py311h06a4308_0
openssl              3.0.15-h5eee18b_0 --> 3.0.16-h5eee18b_0

```

Downloading and Extracting Packages:

```

openssl-3.0.16      | 5.2 MB | : 0% 0/1 [00:00<?, ?it/s]
certifi-2025.1.31   | 163 KB | : 0% 0/1 [00:00<?, ?it/s]

ca-certificates-2025 | 129 KB | : 0% 0/1 [00:00<?, ?it/s]
certifi-2025.1.31    | 163 KB | : 10% 0.0984012204057609/1 [00:00<00:00, 1.05s/it]

openssl-3.0.16       | 5.2 MB | : 1% 0.011930375800648256/1 [00:00<00:09, 9.26s/it]

ca-certificates-2025 | 129 KB | : 100% 1.0/1 [00:00<00:00, 8.74it/s]
certifi-2025.1.31    | 163 KB | : 100% 1.0/1 [00:00<00:00, 1.05s/it]

```

```

Preparing transaction: done
Verifying transaction: done
Executing transaction: done

```

```

!conda --version
!python --version

🔄 conda 24.11.1
   Python 3.11.11

```

```
!conda config --set channel_priority strict
```

```
!conda install -c conda-forge pyfmi --yes # Install the key package
```



```
Preparing transaction: done
Verifying transaction: done
Executing transaction: done
```

✓ BPL_STEM_AIR_Perfusion setup

Now specific installation and the run simulations. Start with connecting to Github. Then upload the two files:

- FMU - BPL_STEM_AIR_Perfusion_linux_om_me.fmu
- Setup-file - BPL_STEM_AIR_Perfusion_explore

```
# Filter out DeprecationWarnings for 'np.float as alias' is needed – wish I could make filter more narrow
import warnings
warnings.filterwarnings("ignore")
```

```
%%bash
git clone https://github.com/janpeter19/BPL_STEM_AIR_Perfusion
```

```
🔄 Cloning into 'BPL_STEM_AIR_Perfusion'...
```

```
%cd BPL_STEM_AIR_Perfusion
```

```
🔄 /content/BPL_STEM_AIR_Perfusion
```

✓ BPL_STEM_AIR_Perfusion - test

Author: Jan Peter Axelsson

Here we show simulations of stem cell cultivation in an aerated hollow fiber reactor. The reactor volume is kept constant and cells recycled, thus the setup is similar to perfusion cultivation.

The model combines rudimentary cell growth and metabolism combined with times series data of of the metabolic rates: q_{Nmax} , q_{Lc} , and q_{O2} , marked with red in the comprehensive plot.

Ref Greuel et al: "Online measurement of oxygen enables continuous noninvasive evalation of human-induced pluripotent stem cell (hiPSC) culture in a perfused 3D hollow-fiber bioreactor", Biotech. Bioeng., 2019.

✓ Setup of the simulation model

```
run -i BPL_STEM_AIR_Perfusion_explore.py
```

```
🔄 Linux – run FMU pre-compiled OpenModelica
```

```
Model for the process has been setup. Key commands:
- par()      - change of parameters and initial values
- init()     - change initial values only
- simu()     - simulate and plot
- newplot()  - make a new plot
- show()     - show plot from previous simulation
- disp()     - display parameters and initial values from the last simulation
- describe() - describe culture, broth, parameters, variables with values/units
```

Note that both `disp()` and `describe()` takes values from the last simulation and the command `process_diagram()` brings up the main configuration

Brief information about a command by `help()`, eg `help(simu)`
Key system information is listed with the command `system_info()`

```
%matplotlib inline
plt.rcParams['figure.figsize'] = [30/2.54, 24/2.54]
```

✓ About the process model

Here a process diagram is shown of the process. Further information about the culutre stored in the model code is extracted.


```
#process_diagram()
```

✓ Simulations

```
# Process parameters to mimic stem cell example
# - study impact of different values maintenatnce qm
```

```
newplot(plotType='Basic')
for value in [1.0e-6, 10.0e-6, 20e-6]:
    par(qm=value)
    simu(350)
```


⚠ Could not find cannot import name 'dopri5' from 'assimulo.lib' (/usr/local/lib/python3.11/site-packages/assimulo)
 Could not find cannot import name 'rodas' from 'assimulo.lib' (/usr/local/lib/python3.11/site-packages/assimulo)
 Could not find cannot import name 'odassl' from 'assimulo.lib' (/usr/local/lib/python3.11/site-packages/assimulo)
 Could not find ODEPACK functions.
 Could not find RADAR5
 Could not find GLIMDA.



```
# Process parameters to mimic stem cell example
# . comprehensive plot including metabolic rates
```

```
par(qm=1e-6)
```

```
newplot(plotType='Comprehensive')
simu(350)
```



```
describe('culture')
```

```
↗ Reactor culture human-induced pluripotent stem cells - hiPSCs
```

```
describe('Vcc')
```

```
↗ Volume of the reactor : 0.04 [ L ]
```

```
describe('N_start')
```

```
↗ Number of cells at start : 50.0 [ 1E6 ]
```

```
describe('N')
```

```
↗ Number of cells : 5223.752 [ 1E6 ]
```

```
describe('G_start')
```

```
↗ Glucose conc at start : 2.5 [ g/L ]
```

```
describe('G_in')
```

↔ Glucose feed conc : 2.5 [g/L]

```
describe('L')
```

↔ Lactate conc : 0.796 [g/L]

```
describe('D0')
```

↔ Dissolved oxygen conc : 6.044 [mg/L]

✓ Appendix

```
system_info()
```



System information

- OS: Linux
- Python: 3.11.11
- Scipy: not installed in the notebook
- PyFMI: 2.16.3
- FMU by: OpenModelica Compiler OpenModelica 1.25.0~dev-51-ge672d09
- FMI: 2.0
- Type: FMUModelME2
- Name: BPL_STEM.Reactor
- Generated: 2024-10-03T11:45:51Z
- MSL: 3.2.3
- Description: BPL - not used
- Interaction: FMU-explore version 1.0.0

Start coding or [generate](#) with AI.